





Clinically  
approved!  
Can be left  
in place for up  
to 14 days

Leukoplast®  
Leukomed® Control  
Clearly controlling the risk of infection

**Leukoplast®**  
Wound care in best hands



**Clear wound visibility<sup>1</sup>**



**Enables early detection of complications**



**Secure bacteria-proof barrier<sup>1</sup>**

**Leukomed Control's** transparency allows visualization of the wound.<sup>1</sup>


It can comfortably and securely stay in place for up to 14 days without the need for dressing removal.<sup>2</sup>

**Transparent, conformable film**

- Enables assessment of visible signs of infection
- Reduces the need for unnecessary dressing changes

**Bacteria-proof film**

- Protects wound from external contamination



**Skin friendly adhesive**

- Wide adhesive edge for high durability
- Secure yet easy to remove

**Absorbent hydropolymer pad**

- Suitable for dry to low exuding surgical wounds
- Pad design aids conformability

Leukomed® Control				
Ref-no.	Size	Pad Size	Items per Unit	Shipping unit
73230-00	5 cm x 7 cm	2 cm x 4 cm	10	1/28
73230-01	7 cm x 10 cm	4 cm x 7 cm	10	1/20
73230-02	8 cm x 15 cm	4.5 cm x 11.5 cm	10	1/20
73230-03	10 cm x 24 cm	6 cm x 20 cm	5	1/20
73230-05	10 cm x 35 cm	6 cm x 31 cm	5	1/10

<sup>1</sup> von Cramon L, et al. A clinical evaluation of a transparent, absorbent, adhesive wound dressing. British Journal of Nursing. 2017 Vol 26 (20): Tissue Viability Supplement

<sup>2</sup> Data on file.

70553-00018-01